

Medlink Cassettes

Hospital Grade Cassettes with Screens

Durable, Tough, Light-tight, Esthetic



Medlink Cassettes Benefits:

- ★ Aviation grade aluminum for outstanding durability of cassette and a patented strong protection layer on screen makes for a robust product to suit any environmental condition.
- ★ Exceptional film-screen contact. Nitrogen imploded open cell P.U. foam system in cassette ensures no deterioration in film-screen contact over extended use.
- ★ Designed to eliminate noise and back-scatter.
- ★ Compatible with all film types.

RARE EARTH - GREEN

Rare earth green screens use terbium activated gadolinium oxysulfide phosphor for increased efficiency. They are designed for outstanding contrast and detail perception and enable reduction in X-Ray dose, with low tube voltage

Screen Type	Speed Class	Benefit
Green 100	100	Outstanding detail.
Green 200	200	Excellent resolution with medium speed.
Green 400	400	Standard rare earth green screen with high intensification and excellent detail.
Green 800	800	Maximum speed for lowest dose.

RARE EARTH - BLUE

Rare earth blue screens use state-of-the-art rare earth phosphor, which filter unwarranted scattered radiation, thus enabling extra sharpness, minimum dosage, and reduction in quantum noise. Phenomenal improvement in X-Ray tube life.

Screen Type	Speed Class	Benefit
Blue 200	200	Incorporates excellent definition with high speed.
Blue 400	400	Standard rare earth blue screen with high intensification and excellent detail.
Blue 800	800	Maximum intensification with lowest dose. Use with half speed X-Ray film will reduce system speed to 400.

Medlink Cassettes

Hospital Grade Cassettes with Screens

Pricing March 2010

Product Code	Description	Dealer Net Price
MK1417C-100BL	MEDLINK 14X17 Cassette W/100 BL Screen	\$109.00
MK1417C-400BL	MEDLINK 14X17 Cassette W/400 BL Screen	\$109.00
MK1417C-400GR	MEDLINK 14X17 Cassette W/400 GRN Screen	\$109.00
MK1417C-800BL	MEDLINK 14X17 Cassette W/800 BL Screen	\$109.00
MK1012C-100BL	MEDLINK 10x12 Cassette W/100 BL Screen	\$70.63
MK1012C-400BL	MEDLINK 10X12 Cassette W/400 BL Screen	\$70.63
MK1012C-400GR	MEDLINK 10x12 Cassette W/400 GRN Screen	\$70.63
MK1012C-800BL	MEDLINK 10x12 Cassette W/800 BL Screen	\$70.63
MK810C-100BL	MEDLINK 8x10 Cassette W/100 BL Screen	\$51.44
MK810C-400BL	MEDLINK 8x10 Cassette W/400 BL Screen	\$51.44
MK810C-400GR	MEDLINK 8x10 Cassette W/400 GRN Screen	\$51.44
MK810C-800BL	MEDLINK 8x10 Cassette W/800 BL Screen	\$51.44
MK2430C-100BL	MEDLINK 24x30 Cassette W/100 BL Screen	\$61.91
MK2430C-400BL	MEDLINK 24x30 Cassette W/400 BL Screen	\$61.91
MK2430C-400GR	MEDLINK 24x30 Cassette W/400 GRN Screen	\$61.91
MK2430C-800BL	MEDLINK 24x30 Cassette W/800 BL Screen	\$61.91
MK1114C-100BL	MEDLINK 11X14 Cassette W/100 BL Screen	\$85.45
MK1114C-400BL	MEDLINK 11X14 Cassette W/400 BL Screen	\$85.45
MK1114C-400GR	MEDLINK 11X14 Cassette W/400 GRN Screen	\$85.45
MK1114C-800BL	MEDLINK 11X14 Cassette W/800 BL Screen	\$85.45
MK1414C-100BL	MEDLINK 14X14 Cassette W/100 BL Screen	\$104.94
MK1414C-400BL	MEDLINK 14X14 Cassette W/400 BL Screen	\$104.94
MK1414C-400GR	MEDLINK 14X14 Cassette W/400 GRN Screen	\$104.94
MK1414C-800BL	MEDLINK 14X14 Cassette W/800 BL Screen	\$104.94
MK1436C-100BL	MEDLINK 14X36 Cassette W/100 BL Screen	\$289.20
MK1436C-400BL	MEDLINK 14X36 Cassette W/400 BL Screen	\$289.20
MK1436C-400GR	MEDLINK 14X36 Cassette W/400 GRN Screen	\$289.20
MK1436C-800BL	MEDLINK 14X36 Cassette W/800 BL Screen	\$289.20
MK612C-100BL	MEDLINK 6x12 Cassette W/100 BL Screen	\$59.21
MK612C-400BL	MEDLINK 6x12 Cassette W/400 BL Screen	\$59.21
MK612C-400GR	MEDLINK 6x12 Cassette W/400 GRN Screen	\$59.21
MK612C-800BL	MEDLINK 6x12 Cassette W/800 BL Screen	\$59.21
MK717C-100BL	MEDLINK 7x17 Cassette W/100 BL Screen	\$68.88
MK717C-400BL	MEDLINK 7x17 Cassette W/400 BL Screen	\$68.88
MK717C-400GR	MEDLINK 7x17 Cassette W/400 GRN Screen	\$68.88
MK717C-800BL	MEDLINK 7x17 Cassette W/800 BL Screen	\$68.88
MK9595C-100BL	MEDLINK 9.5x9.5 Cassette W/100 BL Screen	\$60.16
MK9595C-400BL	MEDLINK 9.5x9.5 Cassette W/400 BL	\$60.16
MK9595C-400GR	MEDLINK 9.5x9.5 Cassette W/400 GRN Screen	\$60.16
MK9595C-800BL	MEDLINK 9.5x9.5 Cassette W/800 BL Screen	\$60.16

Add \$8.00 for Window Cassette

Prices and specifications are subject to change with 30 day notice