

The Portable Solution



JADE Mobile




MEDLINK IMAGING
a viewworks company

The JADE Mobile Solution

JADE is a powerful, innovative, highly versatile portable X-ray system suitable for a variety of clinical uses. With a simple, easy-to-use operator console, two-step foldable stand and auto lock system, JADE is one of the most user-friendly portable X-ray systems on the market today. When combined with the award winning VIVIX-S 1417W panel, an unparalleled pre-packaged mobile digital diagnostic system was created.



Unit: mm(inch)



JADE Unit Features

1417W Panel Features

<ul style="list-style-type: none"> • 4kW High Frequency X-ray Generator 	<ul style="list-style-type: none"> • Model FXRD-1417WB
<ul style="list-style-type: none"> • Maximum 120kVp / 0.1 ~250mAs 	<ul style="list-style-type: none"> • General Radiography Application
<ul style="list-style-type: none"> • 0.001 to 10 sec Exposure Time Range 	<ul style="list-style-type: none"> • High Spatial Resolution with 140µm pixel array
<ul style="list-style-type: none"> • 119 lbs including X-ray unit, collimator and stand 	<ul style="list-style-type: none"> • 14" x 17" (35cm x 43cm) Image Size
<ul style="list-style-type: none"> • Maximum height of 7' 5" 	<ul style="list-style-type: none"> • 2,560 x 3,072 pixel size
<ul style="list-style-type: none"> • Collimator with 30 second LED lamp timer 	<ul style="list-style-type: none"> • 14 bit Gray-scale
<ul style="list-style-type: none"> • Exposure Hand Switch with 15' cord length 	<ul style="list-style-type: none"> • Lithium Ion 4000 mAh battery – 4h (capturing)
<ul style="list-style-type: none"> • Foldable, two-step stand 	<ul style="list-style-type: none"> • Stable and Reliable Automatic Exposure Detection
<ul style="list-style-type: none"> • Vehicle loadable 	<ul style="list-style-type: none"> • 2 sec (wired) / 4.5 sec (wireless) Acquisition Time
<ul style="list-style-type: none"> • Wheel lock 	<ul style="list-style-type: none"> • Panel Dimensions: 15" x 18" x .6"
<ul style="list-style-type: none"> • Automatic tube arm lock at any angle 	<ul style="list-style-type: none"> • Inside AP- Fast Wireless Communication
<ul style="list-style-type: none"> • User Programmable APR, (up to 9 APR settings) 	<ul style="list-style-type: none"> • Water and Dust Resistance Rating
<ul style="list-style-type: none"> • Dual Focal Spot of 0.5/1.5mm 	<ul style="list-style-type: none"> • Dell Latitude Laptop
<ul style="list-style-type: none"> • USB interface & Blue-tooth 	<ul style="list-style-type: none"> • VXvue Viewer Software Program
<ul style="list-style-type: none"> • 2 year Manufacturer Warranty 	<ul style="list-style-type: none"> • 5 year Panel Manufacturer Warranty